



DentalMart[®]
dentalmartproducts.com



DENTAL PRODUCTS



ADHESIVE TIP APPLICATOR

Used to stick small device like ceramic veneer, inlay and orthodontic bracket, special made glue tip with strong stick force but not remains left over, well-designed, convenience to use to avoid getting slipped from the forceps.

Packing : 20pcs/box

MICRO APPLICATOR TIPS

Bendable and disposable applicator for cavity lining, bonding, etching, fluoride varnishing saving medicine and cost. Different sizes of micro applicator for your choice, eco-friendly, fits for medical requirements.

Length : 10cm

Three different sizes in tip:

Regular size (2.0 mm)	Color	:	Pearl Green
Fine size (1.5 mm)	Color	:	Pearl Yellow
Super fine size (1.0 mm)	Color	:	Pearl White

Packing : 100pcs/Pack



AUTOCLAVABLE BIB CLIP

Autoclavable, silicone material for the colorized part, stainless material for the head part.

SPECIFICATION

Length of silicone material	: 397mm
Head part length	: 51mm
Packing	: 1pc./pkt.
Color	: Light Orange, Light Green, Light Yellow, Light Purple, Light Blue

BUR HOLDER BOX

Endo box is a storage unit for instrument used during root treatments. Made of alloy plastic, which is durable and efficient to use. Autoclave up to 135°C and can be autoclaved up to 3000 times.

Autoclavable can be heated to 135°C

Specification	: 107X49X55 mm
Color	: Light grey, dark brown
Function	: High-speed,expansile bur holder box and gutta percha point
Packing	: 1pc/box





Specification : (L) 132X90X17mm
(M) 118X85X17mm
(S) 97X70X15mm
Packing : 2pcs/set
Remark : Frosted Surface

CHEEK RETRACTOR PK/2 AUTOCLAVABLE

Structure, Performance and Applicable Scope of Product

Our product is composed of C-shaped connecting rod and retractors on both sides with good intensity and ductility. Different from traditional cheek retractor, our product is made by special technique. The frosted surface with half-clear white color brings good texture and hand-touching, also it is anti-skid. When oral treatment, put cheek retractor on both corner of patient's mouth, and open the mouth adequately to reduce patient's fatigue. Cheek retractor is a assistant tool for doctor's treatment.

Product Instruction

1. Confirm the cheek retractor is clear before using.
2. Choose appropriate size according to patient's mouth during treatment.
3. Turn connecting rod down, bend it with thumb, middle finger and ring finger.
4. Finish one side through cooperation of both hands, then finish another side. Adjust it then.

Product Application Notice

1. Product contraindications : none
2. The product is prohibited to use when product is in a damaged or deformation condition.
3. The suggestion sterilize degree is 121 degree.

CHEEK RETRACTOR Y TYPE

Latex free, flexible and comfortable Cheek Retractor provides retraction of lips and cheeks during dental treatment. Cheek Retractors are durable and long handle ensures ease of application. The smooth surface makes patient very comfortable. Autoclave up to 134°C

Autoclavable to 134°C

Packing : 2pcs/box
Color : Clear



COMPULE DISPENSER GUN

Compule Dispenser Gun is used for conveying ceramic veneer bonding resin; it can be reused and resistant to 121°C high temperature. The bullets material can avoid the light, and is a safe injection resin material. Bullet is disposable, and can prevent cross infections.

Specification : 150X100mm

COTTON ROLL DISPENSER

PRESS TYPE

Easy to use and refill. Top cover can be opened for putting into cotton rolls. To get the cotton rolls, just push down the top cover lightly.

Specification : 120X55X97mm
Packing : 1pc/box
Color : White, Blue

DRAWER TYPE

Easy to use and refill. Hygienic, stable and anti-skid. To get cotton rolls just pull out the drawer.

Specification : 108X53X102mm
Packing : 1pc/box
Color : White, Blue, Red, Yellow



PRESS TYPE



DRAWER TYPE



DENTAL SECTIONAL CONTOURED MATRICES KIT

The Kit Contains :

1. One piece of medical forceps.
2. 3 pcs. of sticky sticks.
3. 150 pcs. of ordinary matrix bands (L/M/S each size is 50 pcs.)
4. 60 pcs. of subgingival matrix bands (L/M/S/ each size is 20 pcs.)
5. 20 pcs. of resin matrix bands (L/S each size is 10 pcs.)
6. 6 pcs. of clamping rings.
7. 20 pcs. of clearance wedges with holes.
8. 20 pcs. of wooden wedges.
9. 20 pcs. of software clearance wedges with four types.
10. 6 pcs. of silicone matrix clips.

DENTURE BOX STANDARD

Specification : 95X74X38mm

Color : White, Pink, Light Blue, Light Green



DENTURE BOX WITH TRAY

Specification : 92X85X57mm

Color : White, Pink, Light Blue, Light Green

DENTURE BOX WITH TRAY

All kinds of colourful denture boxes can protect the Dentures, impression trays, retainers and keep from the dust, keep clean and healthy.

Specification : 86X78X78mm

Color : Dark Green, Green



DISPOSABLE DENTAL INSTRUMENT KIT 7-IN-1

Ideal for the diagnosis and treatment of dental, avoid of all kinds of diseases because of cross-infection, Sterilized by ethylene oxide strictly, two years valid date, including : mouth mirror, probe, forceps, bibs, tissue, tray, cotton balls.

Packing : mouth mirror, probe, forceps, bibs, facial tissue, tray, cotton balls. (1 each)



FACE SHIELD WITH FRAME

Face shield with eyeglass, protected from fluid sprays, splatters and splashes, package coming with 1 frame and 10pcs clear detachable shields, lengthen nose pad, frame can be used in conjunction with eyesight correction glasses.

Packing : 1 Frame+ 10 Shield/box
Color : Blue, White

FLORIDE FOAM TRAY

Perforated trays are made from strong, rigid material to generate consistent impressions. These trays are distortion-free, pre-shaped, readily trimmed are flame softened for adjustment. Handles are non-slip grips and are compatible with a wide variety of impression materials.

Specification	Size	Color	Packing
193x65mm	Large	Blue	50 trays per/pack
180x60mm	Medium	Yellow	
160x55mm	Small	Green	



IMPRESSION TRAYS - AUTOCLAVABLE

Perforated trays are made from strong, rigid material to generate consistent impressions. These trays are distortion-free, pre-shaped, readily trimmed, are flame softened for adjustment. Handles are non-slip grips and are compatible with a wide variety of impression materials. Available in disposable and Autoclave material up to 121°C


01~Large upper; 02~Large Lower; 03~Medium upper; 04~medium lower;
 05~Small upper; 06~Small lower; 07~upper left,lower right;
 08~upper right,lower left; 09~upper anterior; 10~lower anterior;

Packing : (1-10) 10pcs/pkt.

Color : Green



INTRA ORAL TIPS



Specifications : 24.2 X 4.2mm

Packing : 50pcs./bag



Specifications : 20 X 2.4mm

Packing : 100pcs./bag

MIXING TIPS

Wide range of Mixing tips allow easy choice as per material suitability with dynamic proportion, the collection suits for Extra – Light/Light/Regular-body impression material, Regular/ Heavy-body impression material etc.



Packing : 50pcs./pkt.

Color : Blue



Packing : 50pcs./pkt.

Color : Green



Packing : 50pcs./pkt.

Color : Brown



Packing : 50pcs./pkt.

Color : Yellow




MIXING BOWL

Flexible plastic bowl is for easy mixing and can be used repeatedly after cleaning. The material is suitable for hard and soft mixture.

Specification	Capacity	Color
130X96mm	690ml	Blue
130X86mm	500ml	Green

MINI ENDO BLOCK

Name	Specification	Color
Root Canal Measuring Table	89X26X32mm	Grey



Mouth Props Set/3

Wedge-shaped Latex Free Mouth prop is ideal for opening the upper & lower arch, for convenient operating. The comfortable and flexible plastic is resistant to high temperature. Autoclave up to 134°C.

Specification	Color	Packing
39.5X39.5X20.4mm	Purple	3pcs/Set
34.5X36.5X17mm	Green	
33X28X17mm	Yellow	

Mouth Prop & Tongue Guard

Easily avoiding patients fatigue by using moth prop, make patients enjoy a relaxing treatment; tongue guard block tongue moving so that no injury would happen during the whole treatment.



SAFETY GLASSES

Protective glass with exquisite design protect eyes from being injured. Anti-fog, anti-hard light, anti-ultraviolet rays with smooth surface provide greater comfort. Anti-Scratch eyewear offers excellent protection and opal view. UV Protection Glass Protect 98% UV rays.

Color : Red, White



SILICON DISPENSER GUN

Dispensing gun works perfectly to dispense dental adhesive silicones, and dental impression material for perfect registration of dental impression. This sturdy device is made of high quality plastic and consists of release lever with a plunger, which coordinates simultaneously to load and remove the cartridges.

Specification : 1:1, 1:10

Packing : 1pc./box

GAS TORCH



SPATULA

Flexible plastic product for easy mixing and can be used repeatedly after cleaning. The material is suitable for hard and soft mixture.

Autoclavable can be heated to 121°C

Specification : 195x32mm


Color : Pink, White, Yellow, Blue

WOODEN WEDGES ASSORTED PACK

Eco-friendly wooden wedges, round tips protect gingival from damage, smooth round edge, square stern. Six different sizes fits for different requirements

Color	L (± 0.5mm)	X (± 0.1mm)	Y (± 0.1mm)	H (± 0.5mm)
Orange	11mm	3.5mm	1.5mm	1.3mm
White	14mm	4mm	1.5mm	1.5mm
Green	14mm	2mm	1.5mm	1.5mm
Yellow	16mm	4mm	1.5mm	1.5mm
Blue	14mm	4mm	2mm	1.5mm
Pink	16mm	4.5mm	2mm	1.5mm

Packing : 100pcs/pack



X RAY FILM HOLDER

Angled with both ends working, in order to take both anterior & posterior films. Helps in taking perfect & accurate pictures. X-ray film holder helps to align with the appropriate tooth structure easily and precisely every time. Saves time and eliminates errors; is patient friendly. Made of sturdy precision moulded plastic.

Specification : 126x25x11mm



Imported & Marketed By :

AVCO DENTAL

NEW DELHI

134, Gr. Floor, Arjun Nagar, Bhisham Pitamah Marg, Opp. Defence Colony, A-Block, **New Delhi**-110003. (INDIA)
Ph. :+ 91-11-24643501/2/3/4, 41633030, Mobile : +91-9900770088
E-mail: avcoindia@gmail.com

MUMBAI

202, Rebello Heritage, 2nd Floor, Above Magaveera Bank, Azad Road, Gundavali, Andheri (E) **Mumbai** - 400069, (INDIA)
Telefax : +91-22-26846616, Mobile : +91-9819049405
E-mail : avcoindia2@yahoo.co.in

BANGALORE

No.141, IInd Floor, Ist Main Raod, Dr. Rajkumar Road, Ist Block, Rajaji Nagar, **Bangalore**-560010, (INDIA)
Phone : +91-80-32003018, 23320055, Fax +91-80-41712563, Mobile : +91-9900770088, 9740001031,
E-mail : avcoindia3@gmail.com