



ENDODONTICS

Microsurgical	F1
Retractors	F8
Explorers & Excavators	F9
Spreaders	F10
Pluggers	F12
Rubber Dam Products	F15
DownPak™ Tips & Accessories	F20





ABOU-RASS MICROSURGICAL EXPLORERS

Hu-Friedy's microsurgical explorers come in various working ends to provide the clinician with useful options when treating simple or complex cavities. They have excellent balance and the resiliency practitioners need for careful yet thorough exploration. These explorers, and many of Hu-Friedy's endodontic microsurgical instruments, were designed in conjunction with Dr. Marwan Abou-Rass.



Apical Explorer

MEX1


Used to initiate the tracking groove in the root tip.



1
Abou-Rass Endodontic Explorer

MAR-EX1

All-purpose explorer that will not bend under pressure.



FLAT

3
Abou-Rass Endodontic Explorer

MAR-EX3

All-purpose 30mm probe for measuring direct root length in search of the apex or any direct measurement.

SINGLE-END



2
Abou-Rass Endodontic Explorer

MAR-EX2

Used to enlarge the opening made by MAR-EX1 into the cortical bone to allow drainage.



4
Abou-Rass Endodontic Explorer

MAR-EX4

Penetrates small osseous defect in cortical bone where MAR-EX1 is too large.



5
Abou-Rass Endodontic Explorer

MAR-EX5

45°-angled explorer may be used to check dental defects, such as root caries, root resorption and fracture lines.



6
Abou-Rass Endodontic Explorer

MAR-EX6


90°-angled, sharp-pointed explorer for detecting apical foramen and apical seal quality.



MICROSURGICAL MIRRORS

Double sided mirrors are manufactured from highly polished stainless steel and come in round or rectangular shapes. Different sizes and shapes allow for ideal access and visibility in all clinical situations. Rhodium-coated front surface mirrors offer sharp, scratch and distortion-free visibility.

VERSATILE OPTIONS



VELVART



RHODIUM COATED



DID YOU KNOW?

Handles of the MM1, MM2 and MM3 are marked with one, two or three bands for easy identification.

* Handle sold separately. See page E7. MBW wrench required for assembly.



ABOU-RASS MICROSURGICAL CURETTES

Useful for tissue degranulation or placing filling material during apical treatments.



ABOU-RASS MICROSURGICAL EXCAVATORS

Long-shanked excavators are especially useful for access into difficult and confined areas.






ABOU-RASS MICROSURGICAL PLUGGERS

Designed with one round tip and one oval tip for retrofill condensation. Straight apical pluggers assure straight and direct access to lingually positioned apices. Angled apical pluggers are useful in narrow isthmus preparations, thin roots or grooved apical preparations.

STRAIGHT



ANGULATED





KHAYAT AND VELVART MICROSURGICAL PLUGGERS

The Khayat and Velvart pluggers are useful for deep canal condensation, and can be used together for this purpose. Designed in conjunction with Drs. Bertrand Khayat and Peter Velvart.

KHAYAT MICROSURGICAL FILLING INSTRUMENT

This Khayat instrument has one angled and one straight spatula end; the straight end is excellent for hard to reach areas in complex deep canal apical procedure. Designed in conjunction with Dr. Bertrand Khayat .



Khayat
Apical
Plugger

PLGBK6

Velvart
Apical
Plugger

PLGVEL56



Khayat
Apical Filling
Instrument

PF18ABK6



ABOU-RASS MICROSURGICAL BURNISHERS

Straight ball burnishers directly access lingual or distal locations with ease and visibility. Angled ball burnishers feature a long shank allowing optimal reach and visibility of the working end.

ANGULATED

STRAIGHT



MICROSURGICAL RETRO FILLING INSTRUMENTS

These retro filling instruments are sized to work with ultrasonic preparation tips.





MICROSURGICAL PLUGGERS/BURNISHERS

Designs feature a small bar burnisher with straight, right, left and universal condensers in three different sizes for versatile usage.





ALLAM RETRACTORS

A set of seven ergonomically-designed retractors with a common handle designed to reduce hand fatigue and enhance visibility. Contact ends facilitate placement over significant structures by retracting the entire flap and providing a large contact surface. Developed in conjunction with Dr. Charbel Allam.



Allam Retractor
Palatal

SRTA01



Allam Retractor
Mandibular Left

SRTA02



Allam Retractor
Mandibular
Right

SRTA03



Allam Retractor
Universal

SRTA04



Allam Retractor
Universal

SRTA05



Allam Retractor
Maxillary Right

SRTA06



Allam Retractor
Maxillary Left

SRTA07



Allam Retractor
Handle

TRMGRIFF

Handle sold separately.

ARENS RETRACTORS

These retractors were designed in cooperation with Dr. Donald Arens.



1 Anterior
Arens
Tissue Retractor

TRARENS1




2 Posterior
Arens
Tissue Retractor

TRARENS2



EXPLORERS

Long shanked explorers used to locate orifices.



DG16

DG16/23


EXDG16

EXDG16/23

#31, #6, #7

EXDG16T

Thick, #6




DG16/17

EXDG16/17

#6

EXCAVATORS

Skanks have extended lengths for easy access and dentin removal.



2
Glick

EXCGL2



31L

EXC31L

#41, #6



31W
Oval Spoon

EXC31W



32L

EXC32L

#41, #6




33L

EXC33L

#41, #6

ENDODONTIC ACCESSORIES

Essential accessories for efficient treatment.



Locking pliers have horizontal and vertical grooves for holding points.



1
Mirror Handle

MH6692K

With millimeter scale.



16 1/2 cm / 6 1/2"
Locking Pliers

EPL1



6" Stainless
Steel Ruler

CLR6




SPREADERS

Spreaders are used to compact filling material during lateral condensation.

D₁=Diameter of instrument 1mm from the tip. These measurements (i.e., .30) appear next to the tips of the instruments shown below.

ISO-SIZED SPREADERS




DID YOU KNOW?

ISO Pluggers have an ISO color-coded band on the handle for ease of identification during usage:
#30-Blue, #40-Black, #50-Yellow, #60-Blue




NICKEL TITANIUM ENDODONTIC SPREADERS

Ni-Ti provides extra flexibility and memory; especially useful for complex treatment sites.



HEAT CARRIERS







PLUGGERS


Used to compact filling material during vertical condensation.
Marked at 5mm intervals to assess penetration depth.

D1=Diameter of instrument 1mm from the tip. These measurements (i.e., .30) appear next to the tips of the instruments shown below.

ISO* SIZED PLUGGERS



RETRO FILLING PLUGGERS



DID YOU KNOW?


ISO Pluggers have an ISO color-coded band on the handle for ease of identification during usage:
#30-Blue, #40-Black, #50-Yellow, #60-Blue




PLUGGERS

D1=Diameter of instrument 1mm from the tip. These measurements (i.e., .40) appear next to the tips of the instruments shown below.

ANTERIOR



POSTERIOR



DID YOU KNOW?

Anterior condensers are 7mm longer than posterior condensers.


Pluggers specifically designed for anterior or posterior usage.



SLEIMAN PLUGGERS / HEAT CARRIERS


Tapered working end facilitate increased pressure in the root canal, achieving an optimal 3 dimension fill. Each tip is sharp in order to penetrate gutta percha and for heat distribution. Designed in conjunction with Dr. Philippe Sleiman.

TAPERED WORKING END



	1	2
RCPSL172	0.40mm	0.74mm
RCPSL275	0.50mm	0.80mm
RCPSL376	0.66mm	0.96mm
RCPSL478	0.78mm	1.20mm
RCPSL577	0.90mm	1.45mm

LUKS



Sleiman Root Canal Plugger .40mm - .74mm	Sleiman Root Canal Plugger .50mm - .80mm	Sleiman Root Canal Plugger .66mm - .96mm	Sleiman Root Canal Plugger .78mm - 1.20mm	Sleiman Root Canal Plugger .90mm - 1.45mm
RCPSL172	RCPSL275	RCPSL376	RCPSL478	RCPSL577

DID YOU KNOW?

Sleiman pluggers feature tapered working ends that facilitate an optimal 3 dimensional fill and function as heat carriers.



FORCEPS/PUNCHES/RUBBER DAM ACCESSORIES

Essential accessories for efficient treatment.

RUBBER DAM CLAMP FORCEPS



4
17½cm/7"


RDF4



Rubber Dam
Forceps
17cm/6½"

RDF

RUBBER DAM FRAMES



Rubber Dam
Frame - Child
13cm/5"

RDCF5



Rubber Dam
Frame - Adult
15cm/6"

RDAF6

RUBBER DAM PUNCHES



Rubber Dam
Punch
16cm/6¾"

RDP



Ainsworth
17cm/6¾"

RDPA

POST AND SILVER POINT REMOVAL FORCEPS



Steiglitz
Straight
14cm/5½"

RFSTR



Steiglitz
(Archer) 45°
13½cm/5¼"

RF45



Steiglitz
90°
12cm/4¾"

RF90




Peet
45°
12cm/4¾"

RFP



RUBBER DAM ACCESSORIES


RUBBER DAM TEMPLATES



Universal
Rubber Dam
Template

RDTUNIV

Adult template
shown. Pedodontic
template also
included.




FDI
Rubber Dam
Template

RDTFDI

Adult template
shown. Pedodontic
template also
included.

RUBBER DAM CLAMP ORGANIZING BOARD




Rubber Dam
Clamp
Organizing
Board

RDCOB

Used to store, sterilize and
organize rubber dam clamps.

Dimensions 4.62" x 2.4" x .375"
(117mm x 61mm x 10mm)

RUBBER DAM KIT



RUBBER DAM KIT

Includes items listed below (also available individually)


	RDKIT
Rubber Dam Clamp Forceps	RDF
Rubber Dam Punch	RDP
Rubber Dam Frame - Adult	RDAF6
Universal Rubber Dam Template	RDTUNIV
FDI Rubber Dam Template	RDTFDI
Rubber Dam Clamp Organizing Board	RDCOB
#0 Satin Steel® Rubber Dam Clamp	RDCM0
#1 Satin Steel® Rubber Dam Clamp	RDCM1
#2 Satin Steel® Rubber Dam Clamp	RDCM2
#7 Satin Steel® Rubber Dam Clamp	RDCM7
#8 Satin Steel® Rubber Dam Clamp	RDCM8
#8A Satin Steel® Rubber Dam Clamp	RDCM8A
#14A Satin Steel® Rubber Dam Clamp	RDCM14A
#212 Satin Steel® Rubber Dam Clamp	RDCM212



RUBBER DAM CLAMPS

A complete line of Satin Steel® Rubber Dam clamps for anterior, premolar and molar applications.


ANTERIOR



00

RDCM00


Flat jaws and high bow for small lower canines and first premolars.



6

RDCM6


Slightly offset jaws for labial caries on anterior teeth.



9

RDCM9


Offset flat jaws and stiff bow for anterior teeth.



9S

RDCM9S

Offset flat serrated jaws for anterior teeth.



212

RDCM212

Flat jaws for labial caries on anterior teeth.




212SA

RDCM212SA

Slightly offset staggered jaws for odd shaped anteriors.


PREMOLARS



0

RDCM0


Flat jaws and high bow for longer premolars.



1

RDCM1


Slightly festooned jaws for upper premolars.



1A

RDCM1A


Deeply festooned jaws for canines and premolars.



2

RDCM2


Flat jaws for premolars.



2A

RDCM2A


Flat jaws for larger premolars.



2AS

RDCM2AS

Flat serrated jaws for larger premolars.



W2

RDCMW2


Flat jaws for premolars.



27N

RDCM27N


Slightly festooned jaws for premolars.



29

RDCM29

Flat jaws for lower premolars.



209


RDCM209

Small, flat jaws with stiff bow for lower molars.




RUBBER DAM CLAMPS


MOLARS - GENERAL




Deeply festooned four-point jaws for molars.




Deeply festooned jaws and extended bow for distal access to molars.




Slightly festooned jaws for molars.




Flat jaws designed for small upper left and lower right molars; paired with RDCM11.




Flat jaws designed for small upper right and lower left molars; paired with RDCM10.




Deeply festooned jaws for partially erupted or irregularly shaped molars.




Deeply festooned jaws for partially erupted or irregularly shaped molars; larger than RDCM14.




Deeply festooned jaws for partially erupted or irregularly shaped molars.




Flat jaws and stiff bow for molars.




Flat jaws for upper right and lower left molars; paired with RDCM25.



Flat jaws for upper left and lower right molars; paired with RDCM24.



Flat jaws for large molars.



Slightly festooned jaws for large molars.



Slightly festooned serrated jaws for large partially erupted or irregularly shaped molars.



CLAMP SET RDSET7

A popular set for basic coverage


Set includes:

- #8 RDCM8
- #7 RDCM7
- #00 RDCM00
- #9 RDCM9
- #2 RDCM2
- #14A RDCM14A
- #8A RDCM8A



RUBBER DAM CLAMPS

UPPER MOLARS




Flat narrow jaws for small upper molars.




Flat narrow jaws for small upper molars.



Slightly festooned jaws for large upper molars.




Deeply festooned jaws for upper molars.



Flat jaws for upper molars.



Offset jaws for upper right molars; paired with RDCM31.



Offset jaws for upper left molars; paired with RDCM30.




Slightly festooned jaws for upper molars.




Slightly festooned jaws for large upper molars.


LOWER MOLARS




Flat narrow jaws for small lower molars.



Flat narrow jaws for small lower molars.



Flat jaws for lower molars.




Slightly reversed festooned jaws for large lower molars.




Flat jaws for lower molars.



Serrated offset jaws for lower right molars; paired with RDCM13A.



Serrated offset jaws for lower left molars; paired with RDCM12A.



Slightly reversed festooned jaws for small lower molars.



Flat jaws for lower molars.

DID YOU KNOW?


The Satin Steel™ matte finish was designed to reduce the reflection of operatory lights, and help reduce eye fatigue.



DOWNPAK™ TIPS AND ACCESSORIES

DownPak tips and accessories are designed for use with the DownPak Cordless Obturation device. Tips are available in various sizes and materials, including Stainless Steel, NiTi and UltraSoft. Tip tapers range in size from .02mm - .11mm, and are autoclavable.

ULTRASOFT TIPS



Tip
Plugger F
ISO 030/.03

DPTFU

Tip
Plugger FM
ISO 040/.05

DPTFMU

Tip
Plugger M
ISO 050/.07

DPTMU

NITI TIPS



Tip
Plugger XF
ISO 030/.02

DPTXFN

Tip
Plugger FM
ISO 040/.02

DPTFMN

Tip
Plugger M
ISO 050/.02

DPTMN



DownPak Cordless
Obturation Device

STAINLESS STEEL TIPS



Tip
Plugger XF
ISO 030/.02

DPTXFR

Tip
Plugger F
ISO 030/.03

DPTFR

Tip
Plugger FM
ISO 040/.05

DPTFMR

Tip
Plugger M
ISO 050/.07

DPTMR

Tip
Plugger L
ISO 050/.09

DPTLR

Tip
Plugger XL
ISO 070/.11

DPTXLR

CUTTING SPOON TIPS



Small
CSS.01

DPTCSS

Large
CSL.49

DPTCSL

ACCESSORIES



DownPak
Instrument
Bender

DPIB



DownPak
Tip
Carrier

DPTC



DownPak
Device
Holder

DPDH



DownPak
Barrier
Sleeves

DPBS



KEY OPINION LEADER KITS

CHIVIAN-ARENS ENDODONTIC SURGERY KIT* **ENDOSXKIT***

One Cassette—Select From:

Large Signature Series [®] 13 Instrument DIN Cassette, Gray	IMDIN131
Large Signature Series 12 Instrument Cassette, Gray	IM4121
23/UNC15 University of North Carolina Expro	XP23/UNC15
1 Explorer	EXD1
5 Scalpel Handle	10-130-05
149 Periosteal	P149
2/4 Molt Surgical Curette	CM2/4
1 Hourigan Retractor	PH1
1 Anterior Arens Tissue Retractor	TRARENS1
2 Posterior Arens Tissue Retractor	TRARENS2
85 Lucas Surgical Curette	CL85
87 Lucas Surgical Curette	CL87
17/18 Indiana University Curette	SIU17/18
13/14 Columbia Curette	SC13/14
18 Iris Scissors Curved, Delicate	S18
Mathieu-Kocher Perma Sharp [®] Needle Holder	NH5074

*Endo SXKIT Cassette Sold Separately.

KHAYAT ENDODONTIC SURGICAL KIT** **IMKHAYAT2**

One Cassette—Select From:

Large Signature Series [®] 18 Instrument DIN Cassette, Blue	IMDIN138S
Large Signature Series 20 Instrument Cassette, Blue	IM4208
Scalpel Handle	10-130-05E
Microsurgical Scalpel Handle, Stainless Steel [®]	SHDPV
Micro Mirror Round, 2.5mm	MM4
NEW! Khayat Apical Plugger with Satin Steel Handle	PLGBK6
NEW! Khayat Instrument with Satin Steel Handle	PFI8ABK6
85 Lucas Surgical Curette, Spoon Shape, Short Satin Steel Handle	CL856
Velvart Condenser with Satin Steel Handle	PLGVL56
NEW! Perma Sharp General Tissue Pliers, 15cm/6"	TP5072
PR-3 Prichard Periosteal with Satin Steel Handle	PPR3S6
13/14 Columbia Curette, Universal with Satin Steel Handle	SC13/146
34/35 Jacquette Scaler, Miniature, Posterior with Satin Steel Handle	SJ34/356
89/92 Carver Discoid-Cleoids, with Satin Steel Handle	CD89/926
DG 16/17 Endodontic Explorer with Satin Steel Handle	EXDG16/176
Cawood-Minnesota Retractor	CRM2
NEW! Perma Sharp Needle Holder, Castroviejo Curved, Round Handle, 14cm	NH5021RC
14 LaGrange Double-Curved Super Cut Scissors	S14SC

**Includes hinged instrument clips.

HAHN ENDODONTIC MICROSURGICAL KIT*

One Cassette—Select From:

Large Signature Series, 13 Instrument DIN Cassette, Blue	IMDIN138
Large Signature Series, 16 Instrument Cassette, Blue	IM4168
Microsurgical Scissors, Straight, Castroviejo, 18mm	SPV
Microsurgical Tissue Pliers, Curved, 18mm	TPSLCOCMN
Microsurgical Needle Holder, Straight, Castroviejo	NHDPV
Microsurgical Scalpel Handle	SHDPV
Microsurgical Periosteal	PH26M
Spahr Periosteal (Optional)	PPSPAHR6
Mini Scalpel Blade	MB69
Labanca Periosteal, #6 Satin Steel [®] Handle	PPLABANCA6

VELVART ENDODONTIC MICROSURGICAL KIT*

One Cassette—Select From:

Large Signature Series 13 Instrument DIN Cassette, Blue	IMDIN138
Large Signature Series 16 Instrument Cassette, Blue	IM4168
Microsurgical Mirror 3mm	MM4
Microsurgical Mirror 5mm	MM5
Microsurgical Mirror 7mm x 2mm	MM7
Microsurgical Mirror 9mm x 3mm	MM8
Microsurgical Needle Holder, 18cm, Straight	NHDPV
Microsurgical Mirror Handle/Micro Scalpel Handle	SHDPV
Mini Scalpel Blade, 12pc/pack	MB67
Mini Scalpel Blade, 12pc/pack	MB69
Micro Surgical Scissors, 18cm	SPV
Micro Tissue Pliers, 18cm, Straight	TPDAPV
Micro Tissue Pliers, 18cm, Angled	TPDSCPV
Micro Tissue Pliers, 18cm, Straight	TPDSSPV

*Add cassette clips and accessories as needed.



KEY OPINION LEADER KITS

MACHTOU ENDODONTIC PRE-TREATMENT KIT*

One Cassette—Select From:

Large Signature Series® 16 Instrument DIN Cassette, Blue	IMDIN168
Large Signature Series 16 Instrument Cassette, Blue	IM4168
IMS® Signature Cassette Clips, 2	I1000
IMS Bur Cushion, Short Lid, Holds 12	IMS-1372S
Articulation Paper Pliers	APF2E
Hollenbach Carver, #6 Satin Steel® Handle	CVHL3S6
Universal Pliers	DPU17
Endodontic Excavator, #6 Satin Steel Handle	EXC32L6
Explorer, #30 Handle, 12.5mm	EXS6XL
Kelly Hemostat, 14mm	H1
5 Front Surface Mouth Mirror	MIR5
#6 Satin Steel European Mirror Handle	MHE6
Qulix® UNC Color-coded Probe, #6 Satin Steel Handle	PCPUNC156
Loesche 2 Composite/Filling, #6 Satin Steel Handle	PFIGML26
Rubber Dam Frame, Adult	RDAF6
#1 Rubber Dam Clamp	RDCM1
#212 Rubber Dam Clamp	RDCM212
#26N Rubber Dam Clamp	RDCM26N
#4 Rubber Dam Clamp	RDCM4
#8 Rubber Dam Clamp	RDCM8
#9 Rubber Dam Clamp	RDCM9
#W8A Rubber Dam Clamp	RDCMW8A
Rubber Dam Clamp Forceps	RDF4
Iris Scissors, Straight, Delicate, 11.5cm	S17
204S Sickle Scaler, #6 Satin Steel Handle	S204SD6
Anesthetic Syringe, Type CW	SYRCWE

MACHTOU ENDODONTIC TREATMENT KIT*

One Cassette—Select From:

Large Signature Series 16 Instrument DIN Cassette, Blue	IMDIN168
Large Signature Series 16 Instrument Cassette, Blue	IM4168
IMS Signature Series Large Parts Box	M-1273
Cement Spatula	CSA6
1 Perry Dressing Pliers	DP1
DG16 Explorer, #7 Satin Steel Colours Handle	EXDG167
5 Front Surface Mouth Mirror, Simple Stem	M5C
#6 Satin Steel European Mirror Handle	MHE6
10 Root Canal Plugger	RCP10
11 Root Canal Plugger	RCP11
8 Root Canal Plugger	RCP8
9 Root Canal Plugger	RCP9
1/2 Machtou Root Canal Plugger	RCPM1/2
3/4 Machtou Root Canal Plugger	RCPM3/4
11T Spreader	RCSD11T
Iris Scissors, Straight, Delicate	S17

VON ARX APICAL MICROSURGERY INSTRUMENT KIT*

One Cassette—Select From:

Large Signature Series® 18 Instrument DIN Cassette, Blue	IMDIN138S
Large Signature Series 20 Instrument Cassette, Blue	IM4208
Signature Series Cassette Clips, 2	I1000
Scalpel Handle	10-130-05E
#15C Scalpel Blades	10-256-15
84 Lucas Surgical Curette, #6 Satin Steel® Handle	CL846
86 Lucas Surgical Curette, #6 Satin Steel Handle	CL866
32L Excavator, #6 Satin Steel Handle	EXC32L6
E1 Excavator (12-9-15), #6 Satin Steel Handle	EXCE16
Halsted-Mosquito Hemostat, Straight, 12cm	H4
Abou-Rass Microsurgical Excavator	MAR-C1
Apical Burnisher/Plugger, Large, Left	MCBLL
Apical Burnisher/Plugger, Medium, Left	MCBLM
Apical Burnisher/Plugger, Small, Left	MCBLS
Apical Burnisher/Plugger, Large, Right	MCBRL
Apical Burnisher, Medium, Right	MCBRM
Apical Burnisher, Small, Right	MCBRS
Microsurgical Explorer	MEX16
Markel Adult Mouth Prop	MP54
Perma Sharp® Needle Holder Castroviejo, Straight	NH5024
Perma Sharp Needle Holder Crile-Wood	NH5038
Colorvue® Probe Tips, 12 pcs (3-6-9-12)	PCV12PT
Colorvue Satin Steel Handle	PH6
Buser Periosteal, #6 Satin Steel Handle	PPBUSER6
3/9A Prichard Periosteal, #6 Satin Steel Handle	PPR3/9A6
Microsurgical Scissors, Curved, Titanium	STMBH
8A XTS® Satin Steel Composite Filling Instrument	TNPF18A
Tissue Pliers Micro, Straight, Titanium	TPASTMBH