


HIGH - TECH INNOVATIVE OIL-FREE DENTAL COMPRESSORS




"OPTIMIZED FOR PERFORMANCE. DESIGNED FOR LIFE"

AEROVA COMPRESSOR

Dome compressors combines **international technology with Indian manufacturing excellence**, giving you **high efficiency**, reliability, and unmatched durability. Dome also offers a complete range of **advanced compressors designed to meet diverse medical and dental needs**,

AEROVA 40



Compressor Head


- "Engineered for pure air,
- Consistent airflow that protects your handpieces.
- Delivers medical-grade performance."




Easy Drain

- "Designed with a quick-drain system.
- Remove moisture effortlessly.
- keeping your tank clean and rust-free."

Bacteria are a constant threat in dental environments, But **Dome compressors** provide an **added layer of protection**. Featuring an innovative **antibacterial coating inside the tank infused with silver nano particles**. This technology **effectively blocks bacterial growth and prevents corrosion**. It guarantees a **hygienic and reliable air delivery**.



Silver Ion Technology



©2009 Agion Technologies, Inc.

AEROLYTE COMPRESSORS

Dome compressors built as per ASTM standards with SS 304 stainless steel, the tank ensures zero rust formation, higher hygiene, and a longer lifespan.

- It delivers pure, contamination-free air that protects delicate dental handpieces and enhances treatment quality.

AEROLYTE 40




1



1




Compressor Head



- "Engineered for pure air,
- Consistent airflow that protects your handpieces
- Delivers medical-grade performance."

Easy Drain



- Designed with a quick-drain system.
- Remove moisture effortlessly.
- keeping your tank clean and rust-free.

- With an ultra-quiet, oil-free compressor head, it maintains a calm clinic environment while providing clean airflow to safeguard both instruments and patients.
- The consistent pressure output supports smooth operation of scalers, syringes, and high-speed handpieces, while the hygienic internal surface prevents bacterial buildup—making it an ideal choice for modern dental practices.

AEROLYTE COMPRESSORS

Dome compressor combines a robust 2.2 HP dual-head design with pure, oil-free air delivery, giving dentists high airflow, low noise, and consistent pressure for superior clinical performance.

AEROLYTE 60



Semi-Autodrain Air Filters

Since cleaner air is supplied to the **Handpiece**, the overall life of dental instruments increases, reducing maintenance and replacement costs.

- Built to medical-grade standards.
- It protects handpieces, maintains steady pressure, and offers the reliable efficiency every dentist expects.

Model	 kw HP	 V Hz	 Pressure bar	 Litres	 db (A)	 Flow Rate L/min	 Size mm L x W x H	 kg
Aerova 40	0.8 / 1.12	220V 50Hz	5 - 8	40	62	135	470 x 380 x 700	45
Aerolyte 40	0.8 / 1.12	220V 50Hz	5 - 8	40	62	135	470 x 380 x 700	45
Aerolyte 60	1.68 / 2.2	220V 50Hz	5 - 8	60	66	135	650 x 430 x 700	53
Aerolyte 50i	1.2 / 1.6	220V 60Hz	5 - 8	50	65	135	580 x 450 x 800	48



DOME HEALTHCARE PVT LTD

No C - 3, 3rd Floor, KSSIDC Industrial Estate,
Yelahanka, **Bengaluru 560064**

TEL - 080 - 4173 6804

MOB - +91 92880 06156

www.dome-healthcare.com