




*Aerolyte Dental
Compressor*



"OPTIMIZED FOR PERFORMANCE. DESIGNED FOR LIFE".

Aerolyte 50i Compressor



This is a **50-litre** medical-grade air compressor, designed and manufactured 100% in India in collaboration with **MGF Italy**.



Compared to other compressors available in the market, this model features patented technology that ensures significantly lower noise levels, making it ideal for clinics and medical environments.

The tank is built as per **ASTM standards entirely from SS 304 stainless steel**, which means zero rust formation, enhanced hygiene, and a much longer lifespan.



Since cleaner air is supplied to the **Handpiece**, the overall life of dental instruments increases, reducing maintenance and replacement costs.



The compressor is also equipped with an **HOUR METER**, which accurately tracks running hours, helping schedule timely service and ensuring consistent performance and long operational life.

Model	 kw HP	 V Hz	 Pressure bar	 Litres	 db (A)	 Flow Rate L/min	 Size mm	 kg
Aerolyte 50i	1.2 / 1.6	220V 60Hz	5 - 8	50	65	135	356 x 490 x 800	40



DOME HEALTHCARE PVT LTD

No C - 3, 3rd Floor, KSSIDC Industrial Estate,
Yelahanka, **Bengaluru 560064**

TEL - 080 - 4173 6804

MOB - +91 92880 06156

www.dome-healthcare.com